

STEP 1

Equipment. Adult Coaxial Ventilator circuit and exp valve, pharma mini port HME, Appropriate sized CPAP mask



STEP 2

Set up adult ventilator tubing as per invasive ventilation. Attach flow sensors to appropriate ports.



STEP 3

Choose adult/ped on Hamilton T1. Ensure pre-op checks are completed. This will include a flow sensor check as per invasive set up.




STEP 4

Attach HME to end of coaxial circuit. Attach **non vented** CPAP mask to HME. Set up head straps as shown below.




STEP 5

In 'Modes' select NIV.



STEP 6

To deliver CPAP in this mode, reduce ΔP support to '0' and select the desired starting PEEP and FiO2.



STEP 7

This setup is compatible with the above mask AND both the 'Avatar' full face masks and standard Laerdal mask + white elbow



STEP 8

If nebulisation is required. Insert Aerogen solo on the **PATIENT** side of the HME as shown. Inserting it incorrectly will saturate the HME.

