

Equipment Selection

1. F&P RT330 Optiflow tubing
2. F&P RT024 extension kit
3. 1L Water for irrigation
4. Armstrong pharma HME filter

Ensure flow req is <25L before opening



STEP 1

Place an HME filter on inspiratory port. Slide humidifier chamber onto its base & attach RT024 extension tubing to the HME & connect to the nearest port on the humidifier



Step 2


Remove humidifier cables from the bag on trolley. Connect yellow-end cable to yellow port & blue-end cable to blue port. These ports are located on the right side of the humidifier.



This may be easier to complete from the back of the trolley

Step 3

1. Insert 3-prong probe of yellow cable into three-prong socket at the vent end of tubing. **2.** Click 'v' shaped blue probe into 'v' socket next to this. **3.** Insert this end of the tubing into the humidifier.



Step 4

At the free 'nasal cannula' end of the tubing, insert the pointed blue probe into the round socket next to the nasal cannula port. Spike the bag of water & hang to allow chamber to fill




Nasal prong selection

Select & connect appropriate size Optiflow nasal cannula based on flow required. Aiming for 2L/Kg/min up to a max of 25L on this circuit setup. Children >12Kg need the larger circuit

	PRODUCT SIZE	APPROX WEIGHT (KG)	Max. flow
		2 4 6 8 10 12 14 16 18 20 22	
	Premature		8 L/min
	Neonatal		8 L/min
	Infant		20 L/min
	Pediatric		25 L/min

Hamilton set up

Turn on the Hamilton. No flow or leak test required. Set weight. Mode set up: Select Adult/Ped for all patients. In "Modes", select HiFlow O2. Set desired flow & FiO2. Turn on humidifier to NIV



O2 consumption calculation

Complete O2 calculation once high flow is established for >5mins. **1.** Select 'system' on bottom of main screen. **2.** O₂ consumption shown on left hand side in L/min. Enter this in O₂ calc on website.

